

## Casa Med Bariatric

The Casa Med Bariatric fully profiling bed is available in two styles:-  
120cm wide bed with grey plastic headboards and the option of clamp or side rails or  
120cm wide bed with ultra beech wooden headboards, complete with metal integral side rails.



### Features & Benefits

- Easy to clean
- Easy to install & reinstall
- Fixed frame very robust and stable
- No tools required for assembly
- Transport bracket aids storage and transportation
- Quick release Dewert motors and control system
- Dewert lockable handset with head tilt as standard
- Detachable end actuator cable to control box connector
- Height adjustable

- Electric backrest, legrest and knee brake
- Dual action headrest and knee brake adjustment on handset, to help reduce shear when profiling
- Manual lower leg section levelling device
- 100mm diameter, 4 brake castors
- 4-section mattress support, 2 motor adjustments, separate for transport
- Three Year Warranty (Parts only)
- Ten Year Warranty - Metal steel tubes/mattress supports (Parts only)

### Technical Specifications

- Total Weight **95kg**
- Outer Dimensions (W×L) **132×213cm**
- Maximum Safe Working Load **317kg/50st**
- Maximum Patient Weight **286kg/45st**
- Mattress Support (W×L) **120×213cm**



Maximum User Weight

### Ordering Details

Product Code	Bed	Headboards	Height Adj.	Mattress Support
21990B BE	Casa Med Bariatric	Grey Plastic	40 - 80cm	Metal Steel Tube
21991B BE	Casa Med Bariatric Low	Grey Plastic	24 - 64cm	Metal Steel Tube
21990BUltra BE	Casa Med Bariatric Ultra	Beech	40 - 80cm	Metal Steel Tube
21991BUltra BE	Casa Med Bariatric Ultra Low	Beech	24 - 64cm	Metal Steel Tube