

# ENERGI RANGE

## Powerchairs

*“Attractive electric blue finish with vivid artwork”*



CE



Crash Tested

The Enigma Energi Lightweight Powerchairs have many excellent features normally found on more expensive chairs yet are surprisingly affordable. This powerchair has an electric blue finish with vivid artwork, and provides great comfort with extra features as standard.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

## Energi



- » Maximum speed of 4mph\*.
  - » Maximum range of 15 miles on a full battery charge\*.
  - » Available in 41cm (16"), 46cm (18"), 51cm (20") and 56cm (22") seat widths.
  - » Designed for indoor use and limited outdoor use such as shopping centres.
  - » Rear wheel drive system for optimum manoeuvrability.
  - » Detachable height and width adjustable armrests with removable PU arm pads to suit comfort of user.
  - » Double cross brace for responsive manoeuvrability and added support.
  - » Half folding and removable backrest for easy storage and transportation.
  - » Tension and height adjustable backrest and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences.
  - » De-clutch operation to enable you to safely freewheel the chair where and when required.
  - » Padded upholstery with choice of four seat widths.
  - » 20cm (8") front castors and 32cm (12.5") rear wheels with solid puncture proof tyres provide durability and low maintenance.
  - » Comfortable seat cushion included as standard.
  - » Swing in / swing out detachable footrests for maximum comfort and convenience when transferring.
  - » PG Drives VR2 controller provides an easy and enjoyable driving experience.
  - » Powerful 320W motors.
  - » Tried and tested industry standard electronics.
  - » Detachable batteries and basket for easy storage.
  - » Built in storage bag for battery charger.
  - » On board charging (2 x 34Ah batteries).
  - » Rear anti tip wheels as standard.
  - » Available in three colour options – Metallic Blue, Red or Black.
  - » This powerchair has been crash tested providing added safety when being transported in a vehicle.
  - » Various optional extras are available such as attendant controls, dual controls, kerb climbers and elevating footrests.
  - » Maximum weight capacity 135 - 180kg (21 - 28st).
- \* Subject to conditions of use.

## Energi +



- » Maximum speed of 4mph\*.
  - » Maximum range of 15 miles on a full battery charge\*.
  - » Available in 46cm (18") seat width.
  - » Stylish oval tubing with an Electric Blue finish.
  - » Excellent slide out battery tray for easy charging.
  - » Designed for indoor use and limited outdoor use such as shopping centres.
  - » Rear wheel drive system for optimum manoeuvrability.
  - » Detachable height and width adjustable armrests with removable PU arm pads to suit comfort of user.
  - » Double cross brace for responsive manoeuvrability and added support.
  - » Half folding and removable backrest for easy storage and transportation.
  - » Tension and height adjustable backrest and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences.
  - » De-clutch operation to enable you to safely freewheel the chair where and when required.
  - » Padded upholstery for additional comfort and support.
  - » 20cm (8") front castors and 32cm (12.5") rear wheels with solid puncture proof tyres provide durability and low maintenance.
  - » Comfortable seat cushion included as standard.
  - » Swing in / swing out detachable footrests for maximum comfort and convenience when transferring.
  - » PG Drives VR2 controller provides an easy and enjoyable driving experience.
  - » Powerful 320W motors.
  - » Tried and tested industry standard electronics.
  - » Detachable batteries and basket for easy storage.
  - » Built in storage bag for battery charger.
  - » On board charging (2 x 34Ah batteries).
  - » Rear anti tip wheels as standard.
  - » This powerchair has been crash tested providing added safety when being transported in a vehicle.
  - » Various optional extras are available such as attendant controls, dual controls, kerb climbers and elevating footrests.
  - » Maximum weight capacity 135kg (21st).
- \* Subject to conditions of use.

### Colour Options (Energi Only)



Metallic Blue    Red (16", 18" & 20")    Black

Leaflet Code: LL202 If you would like to view this product please contact your local stockist, details below.

## ENERGI RANGE

CODE: PCA16 / PCA18 / PCA20 / PCA22 (Energi) PCPI8 (Energi +)

Max Weight Capacity: 135/180kg (21/28st)



Length:  
118cm (46") - E  
114cm (45") - E+



Height:  
94cm (37") - E  
96cm (38") - E+



Width:  
61/66/71/76cm  
(24/26/28/30") - E  
66cm (26") - E+



Seat Width:  
41/46/51/56cm  
(16/18/20/22") - E  
46cm (18") - E+



Seat Depth:  
40cm (16") - E  
40cm (16") - E+



Weight:  
56 - 58kg (123 - 128lb) - E  
60kg (132lb) - E+