

# GEMINO

A new state-of-the-art rollator



- Easy to operate
- Self-locking folding system (Patent pending)
- Unique height adjustment (Patent pending)
- Unique stability
- Excellent design
- Easy to transport
- Low weight

# Gemino

## Easy to adjust



Correct height adjustment of the rollator is important for achieving the best upright walking position. Height adjustment is performed without use of any tools; pull out the handle, make the adjustment, release the handle. The push handle height is automatically and safely locked in the chosen position (Patent pending).

## Easy to transport

Low weight, only 7.1 kg, and small external dimensions when folded make Gemino easy to transport. Gemino has a self-locking folding mechanism that works in both folded and unfolded position (Patent pending).

## Easy to operate



With its stable and firm construction, Gemino increases mobility and encourages activity.

The seat sling is a mesh design, making it easy to see when you walk.

Anatomically designed handgrips with brakes make the rollator and brakes easy to operate.

## Easy maintenance and reuse

All surfaces are smooth and therefore easy to clean. The seat sling is made of plastic and dries quickly after cleaning.

Full range of functional accessories.

## Environmentally friendly

Gemino is designed for a long product lifetime, which is one of the best ways to avoid waste and spare the environment. The main constituent of the rollator is aluminium with a highly durable hard anodising surface. All the materials and parts are prepared for easy and safe dismantling and sorting for recycling.

## Specification

Specification	Dimensions Gemino 30	Dimensions Gemino 30 S
Total height, from surface to top of push handles	78-100 cm	65-77 cm
Total length	65 cm	68 cm
Total width	60 cm	55 cm
Width between handles	47 cm	42 cm
Seat height	62 cm	48 cm
Height folded	80 cm	69 cm
Length folded	65 cm	68 cm
Width folded	23 cm	23 cm
Turning diameter	84 cm	79 cm
Front/rear wheels	20 x 3.6 cm	20 x 3.6 cm
Total weight	7.1 kg	6.6 kg
Maximum user weight	150 kg	125 kg
Maximum load in basket	5 kg	5 kg
Recommended user height	150-200 cm	125-165 cm



## Excellent design

The things we surround ourselves with are important for how we are perceived. Gemino's attractive appearance appeals to all user groups.

Handicare's design work not only emphasises giving the product an attractive appearance, but achieving good, self-explanatory and practical features. Gemino's many excellent design details include built-in reflectors that make the rollator visible in the dark from every angle. Another feature is the basket, which is standard on Gemino and can be detached and used as a shopping bag, complete with functional pockets.

## Gemino 30 and Gemino 30 S

Gemino 30 is recommended for people from 150 cm to 200 cm tall, weighing up to 150 kg. Gemino 30 S is narrower and recommended for people from 125 cm to 165 cm tall, weighing up to 125 kg.

## Product and accessories

Article number	Description
7100000	Gemino 30
7100100	Gemino 30 S
7104006	Basket, advanced, Gemino 30
7109060	Umbrella, small, black
7109010	Crutch holder
7109040	Tray, Gemino 30
7109041	Tray, Gemino 30 S
7109020	Back support, Gemino 30
7109021	Back support, Gemino 30 S