

 **QUICKIE**<sup>®</sup>  
Live without limits.

# **JIVE<sup>®</sup> F XL**



# Expanding the Jive Family with the front wheel drive XL version

## Accessibility

With a minimum seat to floor height of only 43cm, access into vehicles and under tables becomes very easy.

## Transportable

The Jive F XL can easily be transferred into vehicles via regular ramps.



**Jay Comfort seat and back** as standard, available in three widths: 56, 60 and 64cm. A headrest can be added at no additional cost.



**4-pole motors**  
Heavy duty 4-pole motors for an improved maximum range.  
14" solid drive wheels and 9" castor wheels for easy manoeuvring and great outdoor performance.



**120AR-netcontroller** as standard for increased functionality.  
Full range of special controls/options.



**Fixed centre mount or swing-away legrests**



Seat sizes:	56,60and64cm(22,23.5,25")	Max kerb height:	5cm (2")
Speed:	(6kph) 3.7mph	Seat tilt (mechanical):	0° / 3° / 6°
Batteries:	60Ah/70Ah	Back angle (mechanical):	-5° - + 12°
Max. range:	Approx. (30 km) 18 miles	Lowest seat to floor height:	43cm (17")
Max user weight:	(240kg) 37.7stone	Min. overall chair length:	116cm (46")
Total weight (with batteries):	(365lbs) 165kg	Max. overall width (with 640 mm seat; no lights):	88cm (34.5")
Max. safe slope:	6° / 11%	Overall drive base width:	64cm (25")

Find out more today! Visit [www.SunriseMedical.com](http://www.SunriseMedical.com) for more information, to find your local dealer, or to request a demonstration.



### Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Jive FXL and other Quickie products through our Facebook page:

[www.facebook.com/QuickieWheelchairs](http://www.facebook.com/QuickieWheelchairs)



Sunrise Medical Limited  
Thorns Road  
Brierley Hill  
West Midlands  
DY5 2BQ  
Tel: +44 (0) 845 605 66 88  
Fax: +44 (0) 845 605 66 89  
[www.SunriseMedical.com](http://www.SunriseMedical.com)