

Sunlift 150

The Sunlift 150 is the ideal hoist for many care settings. It enables carers to perform safe and secure client transfers in a wide range of lifting situations, including lifts from the floor.

Safe working load	150 kg
Max total length (with wheels facing outwards)	1,215 mm
Min total length (with wheels facing inwards)	1,155 mm
Max width with legs closed	665 mm
Max height (without load)	1,930 mm
Min height (without load)	1,255 mm
Max spreader bar height (without load)	1,680 mm
Min spreader bar height (without load)	450 mm
Turning radius with legs open	1,380 mm
Legs open:	
Max width (external)	1,205 mm
Max width (internal)	1,120 mm
Legs closed:	
Max width (external)	625 mm
Max width (internal)	545 mm
Total leg height	100 mm
Ground clearance	20 mm
Castor wheels	75 mm
Rear wheels	100 mm
Mast, boom and electric components assembly weight	19.3 kg
Base assembly weight	15.1 kg
Batteries weight	2.9 kg

Coopers Product Code	TCES Code Number	Coopers Product Name
OP-G150E	-	Sunlift 150



Sunlift 175

The Sunlift 175 is an ideal hoist for the NHS and community equipment store. With its high lifting capacity and lightweight design the 175 works well in almost all care settings.

Safe working load	175 kg
Max total length (with wheels facing outwards)	1,340 mm
Min total length (with wheels facing inwards)	1,265 mm
Max width with legs closed	665 mm
Max height (without load)	1,985 mm
Min height (without load)	1,425 mm
Max spreader bar height (without load)	1,755 mm
Min spreader bar height (without load)	550 mm
Turning radius with legs open	1,500 mm
Legs open:	
Max width (external)	1,220 mm
Max width (internal)	1,140 mm
Legs closed:	
Max width (external)	625 mm
Max width (internal)	545 mm
Total leg height	115 mm
Ground clearance	15 mm
Castor wheels	75 mm
Rear wheels	100 mm
Mast, boom and electric components assembly weight	22.65 kg
Base assembly weight	18.1 kg
Batteries weight	2.9 kg

Coopers Product Code	TCES Code Number	Coopers Product Name
OP-G175	-	Sunlift 175

