

Party on  
Saturday  
at Brian's...

# MY PUMA MOVES WITH ME



I stay  
connected

**PUMA20**



handicare

## It's a brand new day...

Naturally, you expect your power wheelchair to take you comfortably wherever you want to go. And offer you all the essential functions you need for an active daily life. In short, a modern mobility product that fits your own lifestyle. This is fully reflected in Handicare's brand new Puma 20. Robust, practical and affordable, this chair offers simply good driving in every situation. And many optional extras.

Puma 20 is designed to empower you. Front or rear wheel drive, this stylish chair delivers the best driving performance in its class. Technically smart and modular in design, Puma 20 can also be easily customised to suit your specific needs. Combine the Sedeo Lite seat frame with your choice of Sedeo seating and electrical options for enhanced comfort and safety. And enjoy a smooth and carefree ride. Puma 20. Moves with you.

...let's get moving

## Simply good driving...

Designed for intensive indoor and outdoor use, Puma 20 is a practical and robust power wheelchair that offers surprisingly good driving characteristics. Front or rear wheel drive, it is easy to drive and control with your choice of renowned electronics. Enjoy a smooth and stable ride at all times thanks to good suspension – with rubber shock absorbers - and comfort seating. You can choose a higher top speed – a respectable 10 km/h – and larger range with the 74 Ah battery option. With strong motors and high torque, this chair is also very manoeuvrable, even at low speed. Together with its compact design and low seat height this makes Puma 20 also very convenient for use in and around the house.

...maximum independence  
for an active daily life



## All the essentials you need...

If you are looking for simply good driving, Puma 20 is the perfect power wheelchair for you. It provides all the essential functionality you would expect from a good mobility product. With the option of much more... Truly modular in design, Puma 20 offers you the flexibility to easily 'upgrade' your chair. Enhance your driving experience by adding high comfort seating options and advanced electrical options. Your Puma 20 comes with the robust Sedeo Lite seat frame. Combine it with your choice of high-quality Sedeo seat cushion, backrest, armrests and legrests. And choose electrical 25° seat tilt adjustment and lights for enhanced seating comfort and safety, also when driving at night.

...and more



## Distinguishing Handicare design...

Handicare has a reputation for quality and safety. Do not expect anything less from your Puma 20. Our uncompromising approach to performance, comfort and ease of use has resulted in a technically smart, modular design with lots of possibilities. Thanks to the robust base and good suspension, Puma 20 offers the best driving experience in its class. While creating a highly functional and durable wheelchair, however, we have certainly not overlooked the importance of aesthetics. Enjoy comfortable and safe driving in style with Puma 20.

... pleasant and carefree driving



## Specifications

DIMENSIONS			FWD	FWD	RWD	RWD
Maximum user weight	kg	stone	136	21.4	136	21.4
Batteries (1)	Ah	Ah	40 / 50 / 60 / 74	40 / 50 / 60 / 74	40 / 50 / 60 / 74	40 / 50 / 60 / 74
Total weight excl. batteries	kg	stone	77.6	12.2	77.6	12.2
Total weight incl. batteries	kg	stone	106 (40 Ah)	16.7	106 (40 Ah)	16.7
Total length incl. legrests	mm	inch	1160	45.7	1095	43.1
Total width 13" drive wheels	mm	inch	610	24	610	24
Total width 14" drive wheels	mm	inch	655	25.8	655	25.8
Seat angle	°	°	0 - 6	0 - 6	0 - 6	0 - 6
Seat depth (effective)	mm	inch	420 - 500	16.5 - 19.7	420 - 500	16.5 - 19.7
Seat width (effective)	mm	inch	420 - 510	16.5 - 20.1	420 - 510	16.5 - 20.1
Seat height excl. cushion (2)	mm	inch	405	15.9	405	15.9
Seat height incl. cushion (2)	mm	inch	455	17.9	455	17.9
Backrest angle	°	°	92 / 97 / 102 / 107	92 / 97 / 102 / 107	92 / 97 / 102 / 107	92 / 97 / 102 / 107
Backrest height	mm	inch	480	18.9	480	18.9
Diameter castor wheels	mm	inch	230 / 255	9 / 10	230 / 255	9 / 10
Diameter drive wheels	mm	inch	330 / 350	13 / 14	330 / 350	13 / 14
<b>PERFORMANCE</b>						
Maximum speed	km/h	mph	6 / 10	3.7 / 6.2	6 / 10	3.7 / 6.2
Range* (ISO 7176-4), 6 km/h	km	miles	34 (74 Ah)	21.1 (74 Ah)	34 (74 Ah)	21.1 (74 Ah)
Range* (ISO 7176-4), 10 km/h	km	miles	29 (74 Ah)	18 (74 Ah)	29 (74 Ah)	18 (74 Ah)
Obstacle height (outdoor)	mm	inch	60	2.4	80	3.1
Ground clearance	mm	inch	60	2.4	60	2.4
Reversing distance	mm	inch	1200	47.2	1700	66.9
<b>ELECTRICAL OPTIONS</b>						
Electrical tilt adjustment	°	°	0 - 25	0 - 25	0 - 25	0 - 25

ISO 7176-19



1. Affects seat height (min. seat height + 25 mm)  
2. Seat height incl. load, 0° seat angle and 13" drive wheels

EN12182

EN12183

\*Range is a theoretical value; the actual range depends on many factors such as user weight, driving programme, use of electrical adjustments and lights, weather and road conditions, condition of batteries and tyres, etc.



## Why Handicare?

We are driven by a passion to make everyday life easier. For disabled and elderly people, but also for those caring for them. It's been that way with us from the start in 1986. That's when Tom Myklebust and two other wheelchair users, all suffering from a paralysing injury, set up 'Rullestolekspertene' (the Wheelchair experts), after their chance meeting at a Norwegian rehabilitation centre. Their vision was to develop and deliver high-quality custom-made wheelchairs for disabled people. To increase their independence and quality of life, and truly make a difference. Today Handicare, as the company is called since 1995, is one of the leading healthcare companies in Europe and Tom continues to be an inspiration to all of us working with him. And a living example of what Handicare is all about. Making people's everyday life easier by identifying with their individual needs and offering a wide range of well-designed, highly functional mobility products, transfer aids, bathroom safety solutions, stairlifts and car adaptations. Why? Because it is our belief that in life every day matters.

*"I believe that making everyday life easier can really make a difference in people's lives. Sharing this passion with so many enthusiastic and dedicated people at Handicare is a privilege and a pleasure."*

Tom Myklebust





# MY PUMA MOVES WITH ME



handicare

[www.handicare.com](http://www.handicare.com)