



Factory Fitted Accessories:

Innov8 Night Pack	INNOV8/NP
Innov8 Cantilever Side Rails (Standard)	INNOV8/SR
Innov8 Cantilever Side Rails (High)	INNOV8/HSR
X-ray Backrest	INNOV8/XB
Innov8 Bed End	INNOV8/BE/(Colour)
IP66 Ingress Protection	INNOV8/IP66
Innov8 15cm Single Castor	INNOV8/SC150
Innov8 20cm Single Castor	INNOV8/SC200
Innov8 12.5cm Double Castor	INNOV8/DC125
Innov8 15cm Double Castor	INNOV8/DC150

General Accessories:

Innov8 Lifting Pole	INNOV8/LP
Innov8 Transfusion Pole	INNOV8/DP
Cromer Cantilever (Bed Cradle)	1369
Innov8 Balkan Beam Traction	INNOV8/TBB
Traction – Foot End	INNOV8/TFE
Monitor Shelf	INNOV8/MS
Oxygen Cylinder Carrier	INNOV8/OC

Height (cm)	Width (cm)	Length (cm)
41.4 to 81.4	101.6 (including side rails)	217 (237 when fully extended)

Note: Minimum quantities apply for IP66, special castors and bed end infill colours.
For easy selection bed accessory and side rail tables, please see pages 41 and 44-45, respectively.

017.

Beds Hospital Beds / Innov8 2000

Product Code:
INNOV8/2000

Key Features:

- Standard ward bed and base platform (accessories can be added)
- Easy to clean and decontaminate, incorporates BioCote®
- Solid bed panels for maximum infection control and ease of cleaning, incorporating ridges for breathability (IP54 standard ingress protection)
- Electrically operated backrest, height and kneebreak with battery backup
- Trendelenburg, reverse Trendelenburg, auto-contour and auto-regression
- 19.3cm under bed clearance for hoist access
- Simple to use Patient Handset and Nurse Control Panel with handset lockout
- Dual sided CPR levers
- Very easy to manoeuvre and work with, reducing back or handling problems
- Removable head/foot ends with colour coding options for estates management
- Modular and easy to maintain, providing total patient comfort
- UK manufacture and support; 3 year warranty (1st year parts & labour, 2nd and 3rd year parts only)
- Maximum user weight 222kg (35st).
- Maximum safe working load 255kg (40st)