

## Wide Heavy Duty Bariatric Mattress

Our 120cm wide heavy duty bariatric mattresses have been designed for use with our Casa Med Bariatric Beds.



Product Code	DH6BUK120 BE	FFCB BE	FFCB160 BE
Max. User Weight	394kg	294kg	294kg
Risk Level	High	High	High
Waterlow Score	20	20	20
Warranty	4 Years reducing***	4 Years reducing***	4 Years reducing***
Foam			
Foam Type	3 Layer Construction	2 Layer Construction	2 Layer Construction
Thickness	17.8cm	17.8cm	16cm
Size	200×120cm	200×120cm	200×120cm
Cover			
Cover	•	Closure fabric protection to 50,000ppm chlorine	Closure fabric protection to 50,000ppm chlorine
Seam	Sewn	Welded	Welded
Fire Retardency	BS 7177:2008 Medium Hazard	BS 7177:2008 Medium Hazard	BS 7177:2008 Medium Hazard
Vapour Permeable	•	•	•
Waterproof	•	•	•
Weight	20kg	20kg	20kg
	Suitable for use with the following beds:- 21990B BE & 21991B BE	Suitable for use with the following beds:- 21990B BE & 21991B BE	Suitable for use with the following beds:- 21990B BE, 21991B BE, 21990BUltra BE & 21991BUltra BE