

WheelDrive compatibility guide

WheelDrive, the next generation power assist from Sunrise Medical, easily fits on almost all common manual wheelchairs - without any extra parts or changes to the wheelchair ('snap-on principle').

Please contact Sunrise Medical if your wheelchair model is not included in this compatibility guide.

WHEELDRIVE



Wheelchair brand / model
Sunrise Medical
Breezy BasiX ² Plus*
Breezy RubiX ² *
Breezy RubiX ² Plus *
Cirrus 4*
Emineo*
Exigo 10*
Exigo 20*
Exigo 30*
Ibis*
Impera Line*
Quickie Argon ² / Sopur Argon ² *
Quickie Argon Ti
Quickie HeliX ² / Breezy HeliX ² *
Quickie Easy 200
Quickie Easy 300
Quickie HeliX ² Comfort
Quickie Helium / Sopur Helium*
Quickie Life i
Quickie Life / Sopur Easy Life*
Quickie Life R / Sopur Easy Life R*
Quickie Life T
Quickie Life RT
Bischoff & Bischoff
Pyro Light

Wheelchair brand / model
Pyro Start
Revolution R2
S-Eco 2
S-Eco 300
S-Top
Triton
Etac
M100
Invacare
Action 4
Adapt SpinX
Invacare XLT
Küschall Champion
Küschall Competition
Küschall K3
Rea Focus
Top end - Terminator Titanium
Meyra
Eurochair Basic
Eurochair Vario
Excel G5
Primat
Service

Wheelchair brand / model
Life & Mobility
Roxx
O4 Wheelchairs
MultiHopper
OlympicHopper
SwayHopper
Panthera
S2*
S2 Swing*
U2*
U2 Light*
TiLite
Aero 2

* Product combination approved by manufacturer(s)

Please note: The above list is an indication of wheelchair models the WheelDrive wheels can be mounted on. However year of build, Centre of gravity setting and wheel size are of influence. And a precondition is that the receiver of the quick release axle is 12mm or 12.7mm in diameter. Therefore it is essential to try and test it to be 100% sure! Or ask your supplier for advice.