



THE MOST POWERFUL MOBILITY
SCOOTER IN THE WORLD,
FOR USERS UP TO 65
STONE AS WELL AS
THOSE SEEKING
ADVENTURE IN THE
COUNTRYSIDE

ON & OFF ROAD MOBILITY SCOOTER

HORIZON MOBILITY

www.horizonmobility.com

3 PHASE AC MOTOR

The world's first mobility scooter powered by a 3 phase AC Induction Motor providing unbelievable torque for the ultimate in performance and control.



Road Wheels



Rear luggage



Suspension seat



Weather Canopy

MAYAN AC Specifications

Model	Mayan AC
Dimensions	
Length	1475mm / 58in
Width	700mm / 28in
Height (maximum with mirrors extended)	1290mm / 51in
Transport Height (without seat)	550mm / 22in
Weight including 75A/hr Batteries	149kg / 325lbs
Maximum User Weight	65 stone, 412kg / 910lbs
Maximum Speed	
	Mode 1 6km / 4mph,
	Mode 2 13kph / 8mph,
Optional Mode 3 private land	24kph / 15mph
Maximum Obstacle Climb	125mm / 5in
Maximum Obstacle Descend	125mm / 5in
Wheels and Tyres	
Front (Pneumatic)	100 / 65 – 9 (Ø 360mm / 14in) (Tubeless)
Rear (Pneumatic)	100 / 65 – 9 (Ø 360mm / 14in) (Tubeless)
Tyre Pressure	2.0 bar / (30 psi)
Batteries and Charger	
Standard Batteries	2 x 75Ah sealed
Optional Batteries	2 x 100Ah
Standard Charger	24v / 8A

Your Local Dealer:

Horizon Mobility Ltd
10 Neptune Business Centre
Tewkesbury Road
Cheltenham GL51 9FB
Tel 01242 223000 www.horizonmobility.com