

Dynamic Mattress System

The Plus II is a high specification system that is ideally suited for use in a Nursing Home environment or the Community.



Digital Pump

Sidhil's new digital pump gives you the best of both worlds - an easy to use control panel that also offers the infinite control of a digital system.



CPR Unit

Easy to use and rapid CPR unit. Also featuring unique "G" connector incorporating safety catch as an alternative CPR mechanism.

Key Features

- 2 cell alternating dynamic pressure therapy, full mattress replacement system.
- Easy to use digital control system incorporating alternating and static modes with variable pressure settings.
- Automatic return from Static to Alternating mode after 1 hour without manual intervention.
- Unique "Nurse Mode" lifts pressure quickly to highest setting for 20 minutes before returning it to original pressure setting.
- Easy to use Transport mode "blanking plate"
- Visible and audible alarm for low or high pressure, alternating failure and power down.
- Top hinge, positioned to aid profiling, reduces "pinching" of cells in the sacral area.
- Two way stretch, waterproof vapour-permeable top cover
- Revolutionary "cell in cell" structure provides for the highest levels of reliability.
- Ideal for use on profiling beds
- Supplied with nylon carry bag

Product Codes

Plus II Dynamic Mattress System:
DYN/DIG/PLUS/2

Technical Specification

Maximum User weight:
180Kg (28 stone)
Risk Factor:
Medium / High risk
Warranty:
2 year warranty (parts & labour)
Cell design:
21 cells throughout the mattress incorporating static head cells
Operating cycle:
12 minute cycle time
Mattress size:
H: 220mm
L: 2000mm
W: 880mm

Accessories

Plus II Dynamic Replacement Pump:
DYNPUMP/PL/SO/SX

Plus II Dynamic Replacement Cover:
DYN/PLUS/2/COVER

Contact our Customer Service Team for more information
Tel: +44 1422 233 007 / E: sales@sidhil.com / www.sidhil.com