

## SARA STEDY ENHANCED STANDING AID



### A VERSATILE STAND AND TRANSFER AID

*Sara Stedy*®, a stand aid that provides mobility-promoting support and encourages more mobile patients and residents to stand up independently.

With an increased inner width, a new seat shape and a safe working load raised to 182 kg, *Sara Stedy* also meets the standing assistance needs of larger

patients. Enhanced features including a pivoting divided seat and pedal-operated chassis legs ensure improved comfort and support, and make working routines easier for carers.

...with people in mind

[www.ArjoHuntleigh.com](http://www.ArjoHuntleigh.com)

# KEY FEATURES



**1 New seat design**

The innovative pivoting seat improves transfer efficiency and patient stability.



**2 Easy-to-open chassis**

Pedal-operated chassis legs adjust easily for optimum access during transfers.



**3 Built-in seat handgrips**

Handgrips integrated in the seat allow the carer to easily turn the support aid.



**4 Better knee support**

Provides extra stability or firm support when a patient cannot fully control leg movements.



**5 Total-lock castors**

New total-lock castors provide complete stability when the aid is stationary.



**6 Lower chassis**

Smaller castors and redesigned legs mean a lower chassis for better access.



| Product Information |                     |
|---------------------|---------------------|
| Weight              | 29.4 kg (64.82 lbs) |
| Max. safe working   | 182 kg (400 lbs)    |
| Total length        | 905 mm (32 ¾")      |
| Total height        | 1051 mm (41 ½")     |
| Height of chassis   | 100 mm (4")         |

|   |                |
|---|----------------|
| Internal width legs open  | 729 mm (28 ¾") |
| Low friction castors, the rear two with brakes  |                |
| Lifter – Protection class IPX4  |                |
| Approvals   |                |
| The product are in compliance with current legislation and applicable product standards |                |

ArjoHuntleigh is a branch of Arjo Ltd Med. AB.

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2011.

**ArjoHuntleigh UK – Head Office**  
 310-312 Dallow Road, Luton, Bedfordshire, LU1 1TD, UK  
 Sales Phone: +44 (0) 1582 745 700  
 Sales Fax: +44 (0) 1582 745 745  
 Rental 24hr Helpline Lo-call Phone: +44 (0) 8457 342 000

Rental Fax: +44 (0) 1582 745 843  
 Email: sales.admin@arjohuntleigh.com

**ArjoHuntleigh UK – Service and Technical Centre**  
 St Catherine Street, Gloucester, GL1 2SL, UK  
 Phone: +44 (0) 8456 114 114  
 Service Fax: +44 (0) 1452 525 207  
 Email: ukservice@arjohuntleigh.com