


zippie[®]
by QUICKIE



RX kidz

A True Paediatric Chair
that Grows with the Child

Fit for Adventures




Centre of gravity and seat height adjustment via vertical axle plate (option)



Adjustable legrests provide 6" of depth and 9" of height adjustment. No dead spots - infinite adjustment.


Room for Growth The Quickie RX Kidz can grow with the child. It features infinitely adjustable legrests and footplates and an optional 2.5cm or 5cm seat depth growth kit.



Optional 2.5cm or 5cm seat depth growth kit. 90° x 90° seating position always maintained through seat and frame growth kit and adjustable footplates

Outrigged for Stability The Quickie RX Kidz is available with swing-away or outrigged front frame. The latter increases forward stability, allowing the child to gain confidence without the risk of falls when leaning forward or transferring from or to the chair.

Technical specification



Frame colour :	Black, yellow, purple, red, blue, green	Frame versions:	Swing-away or outrigged fixed
Upholstery colour:	black	weight of heaviest part:	4.4kg
Max user weight:	75kg (11.8 stone)	Front wheel diameter:	10, 12, 15cm (4", 5" and 6")
Product weight:	from 13.4kg	Rear wheel diameter:	30, 50, 55 and 60cm (12, 20, 22 and 24")

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



PB_Zippie_RX_Kidz_REV2

Sunrise Medical Limited
 Thorns Road
 Brierley Hill
 West Midlands
 DY5 2LD
 Tel: +44 (0)845 605 66 88
 Fax: +44 (0)845 605 66 89
www.SunriseMedical.co.uk