



# TravelMate



**Are you ready to go?**

**AMIGO MOBILITY**

0800 781 1200

[www.amigoscooters.co.uk](http://www.amigoscooters.co.uk)



# TravelMate

It's portable--and affordable!

0800 781 1200

www.amigoscooters.co.uk

## Features

The Amigo TravelMate is a must-have for the active traveler who is always on the go!

The most compact model in the Amigo line, the front drive TravelMate has a tight turning radius for maneuvering small spaces, and weighs a very manageable 53 lbs (24 kg) once the batteries and basket are removed. The TravelMate also folds easily with a simple press of a foot pedal for storage in trunks or compartments. Let us help you get going today!



Spare battery under seat



Speed range switch:  
High speed range 0-4 mph (0-6.4 kph)  
Low speed range 0-3 mph (0-4.8 kph)

Available colors:



Lime



Pacific



Cinnamon

<b>Model</b>	580003
<b>Travel surfaces</b>	Indoor/outdoor hard flat surfaces
<b>Drive System</b>	Front drive, hub motor
<b>Platform</b>	One piece, aluminum
<b>Weight capacity</b>	250 lbs (113 kg)
<b>Speed</b>	0-4 mph (0-6.4 kph)
<b>Batteries</b>	Standard 24V, (2) 9 AH airline-safe, sealed lithium-ion option
<b>Battery capacity</b>	9.5 miles (15.2 kilometers)
<b>Charger</b>	Universal 2A offboard
<b>Seat</b>	Black vinyl low-back folding
<b>Front tire</b>	8" (22 x 5 cm) solid rubber, flat-free
<b>Rear tires</b>	9" (23 x 7 cm) Trimline solid foam, flat-free
<b>Turning radius</b>	31" (78.7 cm)
<b>Underbody clearance</b>	4" (10.1 cm)
<b>Length</b>	40.5" (102.8 cm)
<b>Rear wheel width</b>	22.4" (56.8 cm)
<b>Height</b>	38" (96.5 cm)
<b>Total weight (without batteries)</b>	53 lbs (24 kg)
<b>Total weight (with batteries)</b>	65 lbs (29 kg)
<b>Warranty</b>	Limited 2-year warranty

## AMIGO MOBILITY

10 Neptune Business Centre  
Tewkesbury Road, Cheltenham, GL51 9FB  
www.amigoscooters.co.uk

Improving Lives Through Mobility® since 1968

