



Taking care of the caretaker

While everybody is thinking of the demands of the patient, who is taking care of the difficulties the caretaker has to overcome each day? Lifting people in and out of chairs and looking after the personal hygiene of disabled people often might cause harm to the attendant and back injuries have become just one of the long-term effects.

The *Aquatec Ocean E-Vip* provides some unusual positioning features for a comfortable working environment without bended backs and heavy loads.



Exceptional wide range of adjustment

No matter whether the caretaker needs the patient at a seating height of more than one meter for easy care or just below 50 cm for effortless transfers, adjusting the height of the *Aquatec Ocean E-Vip* only takes the push of a button. The 35° tilted position which can be electrically adjusted too allows a safe seating position even for patients with lower upper body control.

Water-tight and machine cleanable

All electrical components of the *Ocean E-Vip* are water-tight sealed and the stainless steel frame makes it possible to clean the complete chair in an institutional washing machine.

Features and Options



Chest-belt included

To ensure the necessary safety during the positioning process the *Ocean E-Vip* comes with a standard chest-belt.



Directional lock

One of the rear castors is equipped with a directional lock for better handling.



Headrest included

The comfortable headrest can be adjusted horizontally and vertically.

Technical data



50 - 105 cm



45 cm



73 cm



101 cm



107 cm - 162 cm
(without headrest)

Aquatec
Ocean E-Vip



0° - 35°



39 kg



150 kg

Aquatec
Ocean E-Vip

Colour



Grey



Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2015 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Ocean E-Vip - EU - 08/2015



Yes, you can.®